

Electron microscopy and nanoscience at Argonne: New opportunities through TEAM

Dean Miller

Materials Scientist & Director,

Electron Microscopy Center

Supported by the U. S. Department of Energy, Office of Science - Basic Energy Sciences, under Contract W-31-109-ENG-38.





A U.S. Department of Energy laboratory managed by UChicago Argonne, LLC

Materials Science Division
Electron Microscopy Center

How will "remote" TEAM impact science at ANL? ... and how can remote operation impact TEAM ??

 Teams of scientists engaged in world class science can access unprecedented capabilities of TEAM





 Remotely located experts can help optimize and facilitate operation of TEAM

Electron microscopy and nanoscience at Argonne

The *Electron Microscopy Center* is the focal point for electron microscopy and scattering at Argonne and for a large user community.






The *Center for Nanoscale Materials* is the home for the development and dissemination of the techniques for synthesis, fabrication, characterization and theory of nanomaterials

The EMC and the CNM work closely together to meet the microscopy needs of their users –from characterization of new materials to "heroic" microscopy studies




CNM Research Themes

NanoBio Interfaces




Create bio-inspired materials and processes for energy transduction

X-ray Microscopy




Create images of the nanoworld with hard x-rays

Nanofabrication




Discover new paths for nanostructured materials, including below 10 nm

Theory & Modeling




Towards the 'virtual fab lab'

Nanophotonics



Understand and control optical energy pathways

Electronic & Magnetic Materials




Understand and control charge and spin-based materials for energy and information transport



EMC Resources and Capabilities for Research

The EMC currently supports six major instruments for scientific research

- TEM & SEM



IVEM-Tandem facility
- unique capabilities for studying ion-solid interactions

- High-resolution SEM
- Analytical TEM

Zeiss FIB-SEM
- focused ion beam
(FIB) w/ SEM



FEI Tecnai F20
- state-of-the-art field-emission gun TEM and STEM





FEI Quanta 400F ESEM
- Environmental SEM with capability for in-situ studies



NanoBio Interfaces – functionalized nanoparticles





Figure 3. TEM image of anisotropic TiO_2 nanoparticles functionalized with biotin and coupled together through avidin. The inset shows a high-resolution image of the coupled nanorods with enhanced contrast seen in the region between the two tips.

Functionalized TiO₂ nanorods, from B. Rabatic et al., Adv. Mater. **18**, 1033 (2006)

Tridentate complex of Ti(IV)_{surface} on a TiO₂ nanoparticle, with pyrroloquinoline quinone, from N.M. Dimitrijevic et al., J. Phys. Chem. B **100**, 25392 (2006)



Key Scientific Questions:

- Statics and dynamics of field distributions (electric, magnetic) and charge/spin transport
- Atomic resolution physical, chemical and electronic structure at hybrid materials interfaces

How TEAM can help:

- Reducing/eliminating damage during imaging and analysis - lower keV with high resolution
- Increased analytical sensitivity
- Small probes with high current
- Space for sample environments & for sample rotation to multiple orientations
- Better contrast for low-Z materials




Electronic & Magnetic Materials: Digital synthesis

Novel behavior in heterostructures probed by digital synthesis: e.g. transparent conducting oxides

Key Scientific Questions:

- Atomic resolution of layer structure and defects
- 3-D structure of interfaces
- field distributions (electric, magnetic) and proximity interactions





Digital synthesis of transparent conducting oxides

Reconstructed Phase

How TEAM can help:

- Atomic resolution and high sensitivity
- Reducing delocalization effects
- 3-D structure— avoid the projection problem





Image reconstruction



Scanning Confocal Electron Microscopy (SCEM)




- Many real world problems require excellent resolution (but not atomic scale) but in thicker materials
- Others require depth resolution which is facilitated by confocal operation
- SCEM provides these capabilities and bridges the gap between HREM and STXM

TEM SCEM

Life Science

Materials Science



> 3 micron thick sections



SCEM


How TEAM can help:

- improve spatial & depth resolution (functions of pre/post specimen illumination and collection angles : α & β)
- New operating modes need to be developed to optimize functionality these are NOT found in standard implementations of existing instruments

Depth of Focus/Field (nm)

10¹

0.1





1000


100

Semi-Angle (mR)

Optimizing TEAM through remote access

Remote experts can help optimize performance of TEAM and troubleshoot issues






The interface for TEAM is intended to be simple and easy-to-use for all users, but occasional higher-level intervention may be needed




New Opportunities through TEAM

Aberration correction is leading to new scientific opportunities




Aberrationcorrecting optics and high stability stage are key steps toward atomic scale resolution

Contrast transfer damping envelope for energy-filtered imaging (chemical contrast)




World's first demonstration of chromatic aberration correction for TEM

TEAM C_c-corrector prototype








 $C_c \& C_s$ corrected



World's first energy filtered TEM using chromatic aberration correction


TEAM C_c-corrector prototype



Energy filtered imaging



Towards 3-D atomic scale chemical imaging



TEAM C_c-corrector prototype



Potential impacts of "remote" TEAM

- provide new scientific opportunities for users through TEAM capabilities
 - a benefit in "passive" mode
 - "active" mode allows research teams to work together as usual yet take advantage of new resources
- new capabilities and optimization of TEAM can be transferred to all users via remote control
 - creates a network of "local" experts



