

# AEM DATA SHEET

Please read the policies and procedures for the use of the **AEM** and complete requirements for one or two specific experiments only. Answer all questions.

**Project Duration** \_\_\_\_\_ **Support** for HVEM/ARM \_\_\_\_\_ Short term (<3 mos.) \_\_\_\_\_ Long term (>3 mos.) \_\_\_\_\_

## Material

What is the material \_\_\_\_\_ Approx. chemical composition \_\_\_\_\_

_____ Thin foil	_____ Bulk	_____ Crystalline	_____ Amorphous
_____ Ceramic	_____ Powder	_____ Metallic	_____ Semiconductor cross-sections
_____ Polymer	_____ Biological	_____ Interfaces	_____ Other _____ Unknown

## Preparation

_____ Ion-mill	_____ Electrochemical polishing	_____ Microtome
_____ Powder dispersion	_____ Extraction replicas substrate	_____ Other

## Nature of Specimen

_____ Self-supporting	_____ 3 mm grid	_____ Be (use for low z microanalysis)
_____ Grid supported		_____ Cu (do not use if sample contains Cu)
_____ Carbon Film		_____ M <sub>o</sub> (use for hot stage above 800°C)
_____ Holey Carbon Film		
_____ Glued to grid	_____ Epoxy	_____ Ag paint _____ Other
_____ Substrate for particle suspension other than carbon		
_____ Carbon-coated	_____ Other coatings	
_____ Beam-sensitive	_____ Radioactive	_____ Toxic
_____ Microtomed specimen on what grid material? _____		

## Morphology

\_\_\_\_\_ Precipitates \_\_\_\_\_ Interfaces & grain boundaries \_\_\_\_\_ Bulk \_\_\_\_\_ Other

## Experiment to be performed

_____ EDXS	_____ UTW (Z>6)	_____ HAD (Z>11)
_____ Chemical microanalysis:	_____ Quantitative	_____ Qualitative
	_____ Bulk (>1wt%)	_____ Trace element (<1wt%)
	_____ Standardless	_____ With standard
_____ Chemical mapping	_____ Channelling enhanced microanalysis	_____ Other (specify)
_____ EELS		
_____ Chemical microanalysis:		
_____ Quantitative	_____ Qualitative	_____ Thickness measurement
_____ EXELFS	_____ ELNES	_____ Other (specify)
_____ Diffraction	_____ CBED _____ Rocking beam	_____ SAD _____ uu diff
_____ Imaging		
_____ Digital image acquisition	_____ Other	
_____ Computer		
_____ EDXS data processing	_____ Image processing	
_____ EELS data processing	_____ Other	
_____ Hot Stage (<1300C) with		
_____ EDXS	_____ EELS	_____ Diffraction _____ Video (provide VHS120 tape)
(Please provide copies of all phase diagrams of subject material & Ta along with your target temperature included with this request)		
_____ Development		
_____ Software	_____ Hardware	_____ Other (specify)

## Background (Please rate your level of competence on a scale: 1 = fully conversant; 5 = beginner)

_____ 200CX operation	_____ STEM operation	_____ Diffraction technique
_____ EDXS operation	_____ EELS operation	_____ KEVEX 8000 routine operation
_____ EDXS data processing	_____ EELS data processing	_____ Computing
_____ KEVEX advanced imaging software		

## Check resources needed:

- 1) Is the microscopist comfortable in analytical electron microscopy? ☐ Yes ☐ No
- 2) Is the microscopist to be trained to operate the AEM? ☐ Yes ☐ No
- 3) Will you need the NCEM to provide an operator to run the AEM for you? ☐ Yes ☐ No
- 4) What electron microscopy classes has the microscopist on this proposal completed?

---